



Chronic Hepatitis B (CHB)


KNOW IT BY THE NUMBERS

Start comprehensive biomarker testing, including **quantitative and qualitative HBsAg**, for more informed clinical decisions¹⁻⁴

- Comprehensive biomarker testing may offer more clinical insights^{5,6}
- **HBV DNA and HBeAg** are complementary biomarkers that guide CHB management—HBV DNA is a measure of viral replication that requires serial monitoring, while HBeAg helps define disease phase and differentiates between active and inactive CHB^{3,7,8}
- **HBsAg** is an important marker of active chronic HBV infection²
 - **Qualitative HBsAg^{9,10}**: Detects presence or absence of HBsAg to assist with diagnosis
 - **Quantitative HBsAg¹⁰**: Informs prognosis and clinical decisions, and tracks therapeutic progress

Comprehensive biomarker testing may offer more clinical insights^{5,6}

Real-world monitoring practices may diverge from the guideline recommendations for chronic hepatitis B (CHB)



Insufficient biomarker testing and delayed treatment initiation can increase the risk of disease progression and serious liver-related outcomes.^{7,8,14-16}

CHB leads to cirrhosis, liver failure, or HCC in **15%–40%** of patients.^{14,17}

Persistent viremia due to inadequate monitoring and lack of treatment adjustment contributes to **ongoing liver injury and fibrosis progression**.^{7,15,16}

In 2022, HBV-related cirrhosis or HCC caused an estimated **1.1 million deaths globally**.⁸

Emerging biomarkers may reveal clinical insights that complement AASLD guidelines⁷

AASLD recommendations are based on HBeAg, HBV DNA, and ALT levels⁷

Quantitative HBsAg testing can help guide management of “grey zone” patients.⁷




Figure used with permission. © 2024 Lim YS. Gray zone of hepatitis B virus infection. *Saudi J Gastroenterol.* 2024;30(2):76-82. https://journals.lww.com/sjga/fulltext/2024/30020/gray_zone_of_hepatitis_b_virus_infection.2.aspx

28%–55% of patients with CHB are viremic, but fall into a grey zone without clear guidance on optimal management and treatment. HBeAg-negative patients are in the “grey zone,” in which HBV DNA or ALT levels are borderline between inactive CHB and immune-active, HBeAg-negative CHB.^{7,19}

- Quantitative HBsAg testing is an emerging biomarker gaining traction in international guidelines (EASL, CASL) due to its predictive and prognostic value, but it remains underutilized in US practice²⁰⁻²²

AASLD=American Association for the Study of Liver Diseases; CASL=Canadian Association for the Study of the Liver; EASL=European Association for the Study of the Liver; HBsAg=hepatitis B surface antigen.

In chronic hepatitis B (CHB):

Comprehensive testing of key biomarkers allows for more clarity when characterizing disease phase and severity, assessing patient prognosis, and evaluating treatment response^{7, 20, 21}

Biomarkers of interest

HBV DNA	Quantitative HBsAg	Qualitative HBsAg	HBeAg	Anti-HBe antibodies
				
Full hepatic function panel:			Noninvasive tests of fibrosis:	
ALT, AST, total bilirubin +/- direct bilirubin, alkaline phosphatase, and albumin ¹			Serum APRI or FIB-4, and elastography, FibroScan, ultrasound, or MRI ¹	

For both cirrhotic and noncirrhotic patients, include HCC surveillance with ultrasound and alpha-fetoprotein.^{7,20}

Anti-HBe=anti-hepatitis B e-antigen; APRI=aspartate aminotransferase to platelet ratio index; AST=aspartate aminotransferase; FIB-4=fibrosis-4; MRI=magnetic resonance imaging.

KEY BIOMARKER: Quantitative and qualitative HBsAg testing

HBsAg testing in CHB management

- HBsAg is a viral surface protein and an important marker of active chronic HBV infection²
- It is a measure of the total transcriptional activity of both covalently closed circular DNA (cccDNA) and integrated HBV DNA in liver cells⁹
- HBsAg is believed to promote immune evasion and disease chronicity⁹

HBsAg tests for CHB:

Qualitative HBsAg	Quantitative HBsAg
Detects only presence or absence of HBsAg in serum via standard immunoassays (eg, ELISA, CLIA) ^{9,10}	Measures serum concentration of HBsAg in IU/mL via automated immunoassay ^{3,20}
Clinical utility^{1,7,10}: • Initial diagnosis of HBV and screening for at-risk individuals • Cannot inform disease activity, treatment response, or prognosis, except in rare cases of HBsAg loss	Clinical utility¹⁰: • Quantitative HBsAg testing may help deliver more comprehensive evaluation—informing prognosis, treatment decisions, and tracking progress toward functional cure*

*"Functional" but not sterilizing cure is achievable and should be defined as sustained HBsAg loss in addition to undetectable HBV DNA 6 months post-treatment.²³

CLIA=chemiluminescence immunoassay; ELISA=enzyme-linked immunosorbent assay.

Quantitative HBsAg and its interpretation

qHBsAg	Below reference (<100 IU/mL)	Low ($100\text{--}<1,000$ IU/mL)	Elevated ($\geq 1,000$ IU/mL)
Interpretation	Suggests inactive CHB ^{7,24} • Risk of viral relapse reduced, varying by ethnicity • Increases the specificity of identifying patients with inactive CHB, but reduces sensitivity to 35%	Suggests inactive CHB ⁷ • May predict spontaneous HBsAg clearance in HBeAg-negative patients with a lower viral load • May require less frequent monitoring	Elevated risk of cirrhosis and HCC ¹⁰
Cirrhosis risk ¹⁰	HR 1 (reference)	HR 1.96	HR 3.5
HCC risk ¹⁰	HR 1 (reference)	HR 3.2	HR 5.4

HBsAg loss (below detection threshold of <0.05 IU/mL) serves as a reliable indicator of reduced viral DNA expression, which is associated with improved clinical outcomes.²⁰

Quantitative HBsAg is an important biomarker in assessing clinical response and functional cure.^{7,20}

Clinical utility of quantitative HBsAg

Prognostic Value

- Predicts risk of liver fibrosis and development of HCC⁷

Identify True Inactive HBV Carriers

- Helps guide treatment approach for patients in an indeterminate stage or “grey zone” (HBeAg-negative with borderline levels of HBV or ALT)⁷

Monitor Response

- Predictor of on-treatment virologic control—helps identify potential candidates to withdraw from therapy (in combination with HBV DNA)²⁵

Determine Probability of HBsAg Loss

- Higher probability of spontaneous HBsAg clearance if <100 IU/mL^{7,26}

HR=hazard ratio.

KEY BIOMARKERS: HBV DNA and HBeAg

Real-world monitoring practices may diverge from the guideline recommendations for chronic hepatitis B [CHB].⁷

HBV DNA is a measure of viral replication and guides treatment decisions, but levels can fluctuate⁷

- Serial monitoring is more reliable than cutoff values
- Levels should be interpreted alongside other clinical variables

HBeAg indicates immune activity, helps define disease phase, and differentiates between active* and inactive CHB^{3,7,8}




Figure used with permission. © 2018 Springer Nature. Yuen MF, et al. Hepatitis B virus infection. *Nat Rev Dis Primers*. 2018;4:18035. <https://www.nature.com/articles/nrdp201835>

*Per AASLD guidelines, immune-active chronic hepatitis B is defined by ALT ≥ 2 X ULN or significant histologic disease with HBV DNA $> 2,000$ IU/mL (HBeAg-negative) or $> 20,000$ IU/mL (HBeAg-positive).⁷

ULN=upper limit of normal.

KNOW THE NUMBERS.

Know the complete picture.

Routine testing of all chronic hepatitis B biomarkers, including **quantitative and qualitative HBsAg**, may provide a fuller picture of patients' disease and better inform treatment decisions.^{3,10,28}

- Quantitative and qualitative HBsAg, HBeAg, and HBV DNA tests are commercially available¹⁰



Scan to explore comprehensive testing
on our website

[KnowHepBTesting.com](https://www.knowhepbtesting.com)

References: 1. Screening and testing for hepatitis B virus infection: CDC recommendations – United States, 2023. Centers for Disease Control and Prevention. March 10, 2023. Accessed September 5, 2025. <https://www.cdc.gov/mmwr/volumes/72/rr/rr7201a1.htm> 2. Clinical testing and diagnosis for hepatitis B. Centers for Disease Control and Prevention. January 31, 2025. Accessed September 5, 2025. <https://www.cdc.gov/hepatitis-b/hcp/diagnosis-testing/index.html> 3. Additional hepatitis B blood tests. Hepatitis B Foundation. Accessed September 24, 2025. <https://www.hepb.org/prevention-and-diagnosis/diagnosis/hepatitis-b-blood-tests/> 4. Understanding your test results. Hepatitis B Foundation. Accessed August 8, 2025. <https://www.hepb.org/prevention-and-diagnosis/diagnosis/understanding-your-test-results/#:~:text=Understanding> 5. Wong RJ. *Gastroenterol Rep [Oxf]*. 2025;13:goaf016. 6. Wong G, Lemoine M. *J Hepatol*. 2025;82(5):918-925. 7. Terrault NA, et al. *Hepatology*. 2018;67(4):1560-1599. 8. Hepatitis B. World Health Organization. July 23, 2025. Accessed September 6, 2025. <https://www.who.int/news-room/fact-sheets/detail/hepatitis-b> 9. Lin CL, Kao JH. *Clin Mol Hepatol*. 2016;22(4):423-431. 10. Ghany MG, et al; 2022 AASLD-EASL HBV-HDV Treatment Endpoints Conference Faculty. *J Hepatol*. 2023;79(5):1254-1269. 11. Pham T, et al. *Med Care*. 2023;61(4):247-253. 12. Zhou Y, et al. *J Viral Hepat*. 2022;29(3):189-195. 13. Ma GX, et al. *Healthcare (Basel)*. 2022;10(10):1944. 14. Fattovich G. *J Hepatol*. 2003;39 suppl 1:S50-S58. 15. Zhang Q, et al. *J Clin Transl Hepatol*. 2021;9(6):850-859. 16. Sun Y, et al. *Clin Gastroenterol Hepatol*. 2020;18(11):2582-2591. 17. Lavenchy D, Kane M. Global epidemiology of hepatitis B virus infection. In: Liaw YF, Zoulim F, eds. *Hepatitis B Virus in Human Diseases*. Humana Press; 2016:187-203. 18. Lim YS. *Saudi J Gastroenterol*. 2024;30(2):76-82. 19. You H, et al. *Infect Dis Immun*. 2023;3(4):145-162. 20. European Association for the Study of the Liver. *J Hepatol*. 2025;83(2):502-583. 21. Coffin CS, et al. *Can Liver J*. 2018;1(4):156-217. 22. Mahajan A, et al. *J Viral Hepat*. 2024;31(11):746-759. 23. Cornberg M. *J Hepatol*. 2020;72(3):539-557. 24. Yang HC. *Clin Liver Dis (Hoboken)*. 2024;23(1):e0195. 25. Hirode G, et al; RETRACT-B Study Group. *Gastroenterology*. 2022;162(3):757-771.e4. 26. Tseng T, et al. *Aliment Pharmacol Ther*. 2015;41(10):949-960. 27. Yuen MF, et al. *Nat Rev Dis Primers*. 2018;4:18035. 28. Vachon A, Osiowy C. *Viruses*. 2021;13(6):951.

Trademarks are property of their respective owners.